

# CL3000 Call Logger Software




## Key features:

- Records 25 years of call history
- Customisable reports
- Export facility
- Graphing and average call response report



The Aid Call "Call Logger" software can be installed on your PC enabling you to download, manipulate and analyse call data from the CP3400 and CP5000. The software provides reassurance to the users' and clients' families and also acts as a management tool to log and report whatever happens to the system, including power failure.



A record of events can be downloaded to include such information as number of visits to a client's room, the times at which the carer arrived and left, if the client leaves their bed or even the room itself.

This is a feature which many sites find invaluable in meeting the requirements of those clients whose needs are more demanding.

**0800 052 3616**  
**sales@aidcall.co.uk**

# CL3000 Specifications

<b>Memory</b>	<ul style="list-style-type: none"> <li>• Approximately 50,000 events. Downloadable to local hard drive.</li> </ul>
<b>Integration</b>	<ul style="list-style-type: none"> <li>• Local PC</li> <li>• Network</li> </ul>
<b>Reports</b>	<ul style="list-style-type: none"> <li>• Reports can be generated on zone, dates, call type, ATX, Carer (IRTX) and call duration.</li> </ul>
<b>Hardware Requirements</b>	<ul style="list-style-type: none"> <li>• Windows XP/Vista/7(32bit)</li> <li>• Ethernet Port</li> </ul>



Aid Call, a business unit of Legrand Electric Ltd  
 Unit 4 Brunel Court, Brunel Road, Newton Abbot  
 Devon TQ12 4PB, Tel: **0845 600 3734** Fax: **01626 332 272**  
 Email: [sales@aidcall.co.uk](mailto:sales@aidcall.co.uk) Web: [www.aidcall.co.uk](http://www.aidcall.co.uk)