


TYPICAL EMERGENCY WIRELESS NURSE CALL SYSTEM



ZNC516 15" Touchscreen Display Panel, 1 metre flex with standard plug fitted. Pre-programmed and fitted to stand ready for deployment.



Deployment of the emergency system can be a temporary one or a temp to permanent one depending on the client's needs. Call points can be fitted on walls with adhesive stickers or screwed for a more permanent fixing. Display panels typically come with the table top mount but can also be fitted to wall easily by local technician if required using a spur as below or the plug provided.



Three variations of call point possible:

- ZNC625 for Push for Emergency
- ZNC630 for Pull for Emergency
- ZNC635 for Pull for Cardiac

Each unit is battery operated and pre-programmed ready to be deployed in minutes.

ZNC630 – Call Point with Display and Emergency Pull-Switch




ZNC625 – Call Point with Display



ZNC635 – Call Point with Display and Blue Cardiac Pull-Switch




Below are some recent examples of systems deployed by Aid Call to assist the NHS with business continuity planning.



Two variations of Pear Push (PP) lead are offered:

- E50310 - Standard one button 2 metre PP lead
- E50315 - Standard one button 2 metre PP lead with in-line safety break
- W07495 - The Linen Clip is used to fasten a Pear Push to the bed sheet when it is not stored in the cradle

E50310 – Standard Pear Push Kit



E50315 – Pear Push plus Safety Break Kit

