

Touchsafe Pro Paging System

Key features:

- Robust and waterproof (IP67) pager
- Vibration and selectable tones
- 2 line, 16 character, backlit display
- Integral 'buddy' nurse present feature
- User re-programmable nurse identity
- Exit alarm
- Belt-clip and lanyard options



Aid Call's Touchsafe Pro Intelligent Pager combines the functionality of a pager and a nurse present 'Buddy' in a robust, waterproof portable unit. Paging messages sent from the Aid Call Nurse Call System are displayed on the clear, backlit display and stored in the memory for later review if necessary.

As well as normal pager functions, Touchsafe Pro can be used to trigger an ATX4000 or ATX5000 call point to give a nurse present (Carer) indication on a display panel. Using the Logon Terminal any pager can be allocated to a staff member using their unique Carer ID. In addition the CP5000 shows the Carer ID location and activity. System activity and Carer ID are also shown on the call logger, allowing management to run reports on individual staff and overall performance.

Accidental removal of the Touchsafe pro portable can be prevented by using a door exit beacon or by programming an out-of-range alert. When the Touchsafe Pro unit passes the exit beacon it alerts the user. If the pager leaves site then it will start to beep and can only be stopped by returning it to the charging rack.

Key Advantages:

Improves resident safety

- Delivers calls directly and therefore more quickly to carers
- Carer accountability is created with Buddy ID facility
- Fully IP rated waterproof to 1m depth so can be wiped or fully washed down in support of infection prevention strategy.

Touchsafe Pro Specifications

Quieter and more responsive care home environment

- Eliminates/reduces call bell noise, a significant issue for residents (particularly those with dementia), families, carers and management alike
- Carers eliminate wasted time walking to check nurse call panels and increasing availability for attending to resident needs.

Protect and leverage your investment

- Exit beacon prevents pager loss from being taken off site
- Waterproof and tough pager case reduces maintenance costs
- Increased information for reviews with resident families, CQC and staff performance management



Charging Rack



Logging Terminal

System	• Radio network - Secure, two way radio network • Frequency - 869 MHz Pan European SRD frequency band • Max number of pagers - 999 • Max base stations with slaves - 17 • Max repeaters - 119 • Exit Beacon - Wireless battery powered with 3 year battery life
Portable	• Environmental - Waterproof IP67 • Dimensions and weight - 65mm x 82mm x 16m, 90g • Display - 2 x 16 characters backlit display • Notification - 10 user selectable bleep tones, silent vibrate and LED 'heartbeat' indicator. • Battery life - 72 hours in normal use • Charge time - 1 hour fast charge, nine hour full charge • Charging method - Multi-bay charging Rack • Battery type - Nickel Metal Hydride, 300mAH
Approvals	• ETSI EN300-220-1 v2.3.1 and ETSI ETS300-683

Legrand Assisted Living and Healthcare

Blyth Workspace, Commissioners Quay, Quay Road, Blyth, Northumberland NE24 3AF.

Tel: 01670 352 371 | Email: sales@aidcall.co.uk | Web: www.aidcall.co.uk

Registered Office: Great King Street North, Birmingham, B19 2LF | Registered in England No. 00115834